

## Sample Submission Guidelines

- **The user has to provide the sequencing kit including the flow cell as per the equipment, chemistry and reads/coverage required for their application in the prescribed format (hyperlink).**
- Share the method/protocol followed to generate the libraries.
- Provide the library pool in 1.5ml low absorbent tube with details of molarity of the library pool.
- Sample tube must be clearly labeled with distinct sample identifier.
- Submit the manifest file (if required) electronically.
- Submit the 'sample sheet .csv' file.
- Provide a formatted pen drive/ hard disk for data collection.